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Instructions for use

Sterile, non-absorbable surgical suture material
Needle-suture combination made of silk

1. DESCRIPTION

SABAsilk® is a non-absorbable, sterile, woven surgical suture material consisting of the organic protein fibroin, extracted from the cocoons of silkworms. These extremely fine cocoon threads are cleaned, spun and precision-woven. This suture material can be coated to improve handling properties and to reduce capillarity. Its outstanding features include high flexibility and good tensile strength. SABAsilk® suture material is available undyed in white or dyed black or blue with physiologically safe dye. It is available in different sizes and lengths with swaged surgical needles of various types and diameters. The packaging of SABAsilk® is marked in light blue. SABAsilk® meets the strict requirements of the European Pharmacopoeia and the EC Directive 93/42/EEC on medical devices.

2. AREAS OF APPLICATION

SABAsilk® is used for the adaptation of wound edges and/or ligature

APPLICATION

The suture material to be used must be selected based on the condition of the patient and the size of the wound.

MATERIAL BEHAVIOR

Silk suture material triggers an initial minimally inflammatory tissue response, followed by incremental envelopment in fibrous connective tissue. Although the silk is not absorbable, its progressive degradation of the proteinaceous silk fibers in vivo can lead to a gradual decrease in stability.

3. CONTRAINDICATIONS

SABAsilk® is not intended for use in the heart, the central cardiovascular system or the central nervous system. It should not be used when preservation of the stability of the suture material is required, as this can deteriorate over a longer period in vivo.

4. WARNINGS

SABAsilk® is intended for single use only. The needles may not be resterilized under any circumstances. Re-use can lead to severe and possibly fatal infections in patients. In accordance with statutory provisions, the manufacturer is not liable for this. SABAsilk® may not be used

- if the expiration date is in the past,
- the sterile packaging is damaged or
- if there is doubt as to whether it is intact.

SABAsilk® may be used only by physicians who are familiar with the use of non-absorbable suture material, as the risk of wound dehiscence depends on the site where the sutures are used and on the suture material used. The material behavior in the tissue must be considered thereby (see above). In addition, the surgeon must have sufficient knowledge of the handling required (e.g. knotting techniques). The surgeon is responsible for complications resulting from an incorrect indication, inadequate surgical techniques or an aseptic error. In all cases it is recommended that the surgery be accompanied by corresponding after-care measures. Like any foreign body, longer contact with saline solutions such as those in the urinary and biliary tracts, can result in the formation of concretions.

As with any other suture material, care must be taken with SABAsilk® to ensure that it is not damaged by incorrect handling. Thus, gripping it with a needle holder or a clamp leaves clear marks on the suture. Such damage has a strong influence on the stability of a suture. The needle must also be handled with the appropriate care. The needle should always be grasped at the beginning of the last third. Grasping the needle at its point can impair penetration and result in breakage of the needle. The same applies to grasping the needle at the swaged end. Caution must also be exercised with regard to puncture injuries from needles. Used needles must be disposed of in a suitable puncture-proof container immediately after use.

5. SIDE EFFECTS

Adverse side effects from the use of SABAsilk® can consist of allergic reactions in patients with sensitivities to silk, as an initial inflammatory tissue response or a temporary irritation in the area of the wound. Like all foreign bodies, the suture material can aggravate an existing infection.

6. STERILITY

SABAsilk® needle and suture combinations are sterilized with ethylene oxide. Resterilization of SABApol® is prohibited! The product must not be used if the sterile packaging is damaged or open. Always dispose of opened, unused products.

7. STORAGE CONDITIONS

SABAsilk® must be stored only in normal environmental conditions and protected from moisture and direct exposure to heat. Store only in the original packaging! Do not use this product after the expiration date!

8. LABELING SYMBOLS

	Manufacturer’s name and address
	Dyed, woven and non-absorbable suture material
	Undyed, woven and non-absorbable suture material
	Dyed, woven, coated and non-absorbable suture material
	Undyed, woven, coated and non-absorbable suture material
	Batch designation
	Sterilized with ethylene oxide
	Do not re-use.
	Do not resterilize.
	Comply with the instructions for use.
	Usable until (year-month)
	Store in a dark location away from sunlight.
	Upper temperature limit for storage
	Do not use if the packaging is damaged.

CE 0123 CE mark with code number of the Notified Body

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